



# Injection technique

Artuvetrin® Therapy is administered by subcutaneous (under the skin) injection. Your veterinarian may give you the option of injecting your pet yourself.

Ask your veterinarian to demonstrate the technique to you. The correct injection technique is shown below.

Shake the vial thoroughly before use.

After use, store the vial in the refrigerator (do not freeze).

1. ▶ Draw up the required volume in the syringe.
  - ▶ Point the needle upwards.
  - ▶ Remove any bubbles by tapping the syringe.
  - ▶ Then press the plunger gently until the first drop of liquid appears.
2. ▶ Lift up the skin at the neck and use your index finger to make a small indentation.
3. ▶ Position the needle in the indentation you have made and push it gently through the skin.
  - ▶ Release the skin at the neck.
  - ▶ Press the plunger to empty the syringe slowly and evenly.

